











ABX Diluent

ABX Diluent

ABX Diluent

ABX Diluent



Store at 18 / 25°C



For in vitro diagnostic use



CE mark 98 / 79 / EC.




Consult instructions for use
www.horiba-abx.com



(Manufacturer)
HORIBA ABX
BP 7290 - 34184 Montpellier
cedex 4 - France

20 L


REF 0901020



LOT

(Exp.)

Buffered isotonic solution
for sheathing and diluting leucocytes,
and for the determination
and differentiation of blood cells,
and the measurement of hematocrit
on ABX blood cell counters



FE