

Buffered isotonic solution for the determination and the differentiation of blood cells, and the measurement of hematocrit on ABX blood cell counters.

Reagent	3	1x 3,4L
---------	---	---------

Erythrocyte lysing agent for white blood cell counting, and differentiation, and hemoglobin determination, on ABX blood cell counters.

Reagent	2	1 x 0.3L
---------	---	----------

Enzymatic solution with proteolytic action for the cleaning of ABX blood cell counters.

Reagent	1	1 x 0,5L
---------	---	----------


 Warning, remove the 4 red caps before use.

 If seal is broken, this Minipack should be treated as a biohazardous waste.

GBC 270A-B

ABX Minipack LMG


REF 0602050



LOT



(Exp.)



For storage and transport only.

R1

R2

R3