

ABX Minolyse

ABX Minolyse

ABX Minolyse



Store at 18 / 25° C



For in vitro diagnostic use



CE mark 98 / 79 / EC.



Consult instructions for use
www.horiba-abx.com



(Manufacturer)

HORIBA ABX

BP 7290 - 34184 Montpellier
cedex 4 - France

Erythrocyte lysing reagent for white
blood cell counting and hemoglobin
determination on ABX blood cell counters.



(Exp.)

LOT

1L

REF 0201010



7003147B

ABX Minolyse

ABX Minolyse

ABX Minolyse