







R3

R2

R1



For storage and transport only



ABX Vetpack
Specially designed for ABX Micros ESV 60,
ABC Vet and Scil Vet abc.

ABX Vetpack



REF 0604052



S/

WASTE

A01A00215 revA

Warning, remove the 4 red caps before use.



REAGENT 1

Enzymatic solution with proteolytic action for cleaning HORIBA Medical blood cell counters.

REAGENT 2

Erythrocyte lysing agent for white blood cell counting and differentiation, and hemoglobin determination on HORIBA Medical blood cell counters.

REAGENT 3

Buffered isotonic solution for the determination and the differentiation of blood cells, and the measurement of hematocrit on HORIBA Medical blood cell counters.

4,2 L

5°C 25°C
Store at 5/25°C

IVD

Tel. + 800 67 14 15 16
www.horiba-abx.com/documentation

Made in France

HORIBA ABX SAS
Parc Euromédecine - Rue du Caducée
BP 7290 - 34184 Montpellier Cedex 4 / France