

ABX **Lysebio**

ABX LYSEBIO are used on
HORIBA Medical blood cell counters
to lyse red blood cell and determine
hemoglobin concentration.



LOT

ABX **Lysebio**

0.4L

ABX **Lysebio**



Store at 15 / 30°C

IVD



Tel. + 800 67 14 15 16
www.horiba-abx.com/documentation

Made in France



HORIBA ABX SAS

Parc Euromédecine
Rue du Caducée - BP 7290
34184 Montpellier Cedex 4 / France

REF 0906013



A01A00201 A