








R3

R2

R1




For storage and transport only



Remove the 3 red caps before installing the pack.

Specially designed for scil Vet abc Plus⁺

scil Vet
abc Plus⁺ pack



REF 0604053



WASTE

A01A00150 rev A

Warning, remove the 4 red caps before use.



REAGENT 1

Enzymatic solution with proteolytic action for cleaning hematological blood cell counters.

REAGENT 2

Erythrocyte lysing agent for white blood cell counting and differentiation, and hemoglobin determination on hematological blood cell counters.

REAGENT 3

Buffered isotonic solution for the determination and the differentiation of blood cells, and the measurement of hematocrit on hematological blood cell counters.

4,2 L

5°C 25°C
Store at 5 / 25°C

IVD



Tel. + 800 67 14 15 16



HORIBA ABX SAS

Parc Euromédecine - Rue du Caducée
BP 7290 - 34184 Montpellier Cedex 4 / France