


Store at 5 / 30°C

Enzymatic solution with proteolytic action for the cleaning of abc and ABX Micros ESV 60 blood cell counters.

REAGENT	1	1 x 0,5 L
---------	---	-----------

REAGENT	1	1 x 0,5 L
---------	---	-----------




 If seal is broken, this Vetpack should be treated as a biohazardous waste.

REF 0604052



ABX **Vet**pack

Specially designed for
ABX Micros ESV 60 and Scil Vet abc Plus



R1

R2

R3

For storage and transport only.

